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**ENGLISH** PLEASE READ THESE INSTRUCTIONS CAREFULLY AND SAVE THEM

**Knee orthosis for post-surgery with graduated movement, telescopic version**

**DECLARATION OF CONFORMITY**  
As qualified manufacturer, ORTHOSERVICE AG declares, under its sole responsibility, that this is a class I medical device that has been manufactured and assessed, according to the EU Regulation 2017/745 (MDR). These instructions have been drawn up in application of the aforementioned Regulation. They are intended to ensure an adequate and safe use of the medical device.  
**TRADEMARKS OF MATERIALS**  
Velcro® is a registered trademark of Velcro Industries B.V.  
**PRECAUTIONS**

We recommend that the pressure exerted by the device does not act on parts of the body with wounds, swelling, or welts. It is advisable to not over tighten the device so as to avoid generating areas of excessive local pressure or the compression of underlying nerves and/or blood vessels. If in doubt on how to apply the device, contact a doctor, physiotherapist, or an orthopedic technician. Read the product composition on the internal label carefully. We recommend to not wear the device in the proximity of free flames or strong electromagnetic fields. Do not apply in direct contact with open wounds.

**WARNINGS**  
It is advisable that the device, designed for the specific indications as set out below, be prescribed by a doctor or a physiotherapist and applied by an orthopedic technician, according to the specific needs of the patient. To ensure effectiveness, tolerability, and proper functioning, the application must be carried out with the utmost care. Never alter the adjustment made by the doctor/physiotherapist or orthopedic technician. The responsibilities of the manufacturer shall lapse in case of inappropriate use or adaptation. The orthosis is made and to be used by only one patient. Should the orthosis be inappropriately used, the manufacturer declines all responsibility, as provided by the regulation for medical devices. In hypersensitive individuals, direct contact with the skin may cause redness and irritation. In case of pain, swelling, welts, or any other adverse reaction, contact your doctor immediately, and in case of a particularly severe adverse event, report the incident to the manufacturer and the competent authority of your country. The effectiveness of the orthopedic product is only guaranteed when all the components are properly in use.

**SELECTION/SIZE**

Item	REF. G1046
Size	universal
Color	black

Fits right and left

**MAINTENANCE**

- Do not bleach  No chemical cleaning
- Do not iron  Do not tumble-dry
- Washing instructions: Rods: wash with a sponge soaked in warm water and mild soap; dry with a cloth. Padding: remove the padding and wash them in warm water with mild soap; allow to dry away from heat sources.

Do not dispose of the product or any of its components into the environment.

**INDICATIONS**

- Post-operative treatment that requires control of knee flexion-extension ligament reconstruction, anterior tibial apophysis transposition, fractures of the distal third of the femur, fractures of the proximal third of the tibia, fractures of the patella, tendon damage to the extensor apparatus
- Post trauma phase of knee sprains
- Varus or valgus osteotomies

**CONTRAINDICATIONS**

Currently no known

**CHARACTERISTICS AND MATERIALS**

- Lateral telescopic element made of light aluminum alloy
- Adjustable telescopic length with easy-to-use buttons on valves
- Valves with strap-loops and structural hook-and-loop system
- Non-metallic adjusting joint, light and intuitively adjusted
  - Flexion: from 0° to +120° with intervals of 10°
  - Extension: from -10° to +60° with intervals of 10°
  - Possibility to lock joint from 0° to +60° with intervals of 10°
- Possible of a block of the articulation of the knee
- MTP Pads to limit brace migration and protect leg from direct contact with hinge; lined only externally, cuttable and perfectly adaptable to leg circumference

**ADAPTING**

- Open all straps and temporarily close them upon themselves.
- Separate the padding from the rods with the straps.
- Follow these steps to adjust the length of the telescopic stays for thigh and calf:
  - raise the button provided on each cuff (fig. A1).
  - Slightly raise the apposite length buttone presente su ciascuna valva (fig. A2)
  - Release the button so that it returns inside the corresponding opening on the metal part.
- Apply the padding such that it fits closely to the leg, closing it at the back with the appropriate Velcro® fasteners: first the thigh padding, then the calf padding (fig. B). In order to position the padding at the optimal height, place the orthosis next to the leg. If necessary the padding may be shortened to adapt it to the patient's legs (shorten the padding along the thinner parts in such a way that the padding coloses properly).
- Model the rods (by bending them) to fit them to the shape of the patients' legs (fig. C).
- Adjust the flexion extension as described (points 1-5, hinge adjustment); perform the same operation on both joints.
- Fix the rods to the pads (fig. D): the rods should follow the median line of the leg and the joints should be aligned with the patella so as to rest on the condyles. To align the rods with the median line of the leg, use the lateral malleolus and the large trochanter as reference points.
- Ensure that the Velcro® straps are perfectly aligned with the corresponding paddings.
- Ensure that the two joints are at the same height.
- Velcro® -fasten the straps after inserting them in the appropriate rings (fig. E-F): fasten first the thigh straps beginning with the one closest to the knee; then proceed by fastening the calf straps. Tighten the straps so as to avoid vertical migration of the guide.

**NOTE:** The bending procedure should be performed by hand, with the help of the grips or laying the orthosis on a rounded foot (chair, edge of a bed); do not use your own knee as a lever.

**ADAPTING FOR THE DOCTOR / TECHNICIAN HINGE ADJUSTMENT**

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