# litecross90/91/92/93

REF. 1090 · REF. 1091  
REF. 1092 · REF. 1093

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## ENGLISH

Please read the instructions carefully and thoroughly

### DECLARATION OF CONFORMITY

As the manufacturer, ORTHOSERVICE AG declares with sole responsibility that this product is a Class I medical product and has been manufactured in accordance with the basic requirements of the directives 93/42/EEC and 2007/47/EEC. These operating instructions were created in accordance with the specification of the above-mentioned directives. They aim to ensure the safe and correct use of the medical product.

### TRADEMARKS OF MATERIALS

Sensitive® is a registered trademark of EUROJERSEY; Velcro® is a registered trademark of Velcro Industries B.V.

### SELECTION/SIZE

Item: REF. 1090						
Size	XS	S	M	L	XL	XXL
Circumf. of waist cm	70/80	80/90	90/100	100/110	110/120	120/130
Front height cm	19	19	19	19	19	19
Back height cm	28	28	28	28	28	28

Colour: C - cream / H - blue / S - black

Item: REF. 1091						
Size	XS	S	M	L	XL	XXL
Circumf. of waist cm	70/80	80/90	90/100	100/110	110/120	120/130
Front height cm	23	23	23	23	23	23
Back height cm	33	33	33	33	33	33

Colour: C - cream / H - blue / S - black

Item: REF. 1092						
Size	XS	S	M	L	XL	XXL
Circumf. of waist cm	70/80	80/90	90/100	100/110	110/120	120/130
Front height cm	23	23	23	23	23	23
Back height cm	55	55	55	55	55	55

Colour: cream

Item: REF. 1093							
Size	XS	S	M	L	XL	XXL	XXXL
Circumf. of waist cm	70/80	80/90	90/100	100/110	110/120	120/130	130/140
Front height cm	14	14	14	14	14	14	14
Back height cm	23	23	23	23	23	23	23

Colour: cream

### MAINTENANCE

- Do not bleach  No chemical cleaning
- Do not iron  Do not tumble-dry

Washing instructions:

#### CORSET

- ensure that the Velcro® fasteners are securely closed before washing the corset
- wash the corset separately
- wash by hand or in the washing machine with mild detergent using an ultra-gentle cycle at a maximum of 30°C
- leave to dry away from direct heat

#### PAD

- Handwash with mild detergent at a maximum of 30°C (we recommend using a sponge)
- leave to dry away from direct heat

Do not dispose of the product or any of its components into the environment.

### PRECAUTIONS

It is advised that the pressure exerted by the product does not act on parts of the body suffering from wounds, swelling or lumps. The product should not be tightened too firmly so as to avoid the development of pressure points or compression of the nerves or blood vessels underneath. It is advisable to wear the product over a garment, avoiding direct contact with the skin. If in doubt about how to apply, contact an orthopaedic technician. The device should not be worn in the vicinity of open flames or strong electromagnetic fields.

### WARNINGS

The product may only be used in the event of medical indication or when prescribed by a physician and must be applied by an orthopaedic technician, the competent authority figure for the application as well as for information on the safe use, according to individual needs. To ensure the effectiveness, tolerability, and correct operation, it is necessary that the device is put on with the utmost care. Any changes to the structure or adjustment of the device must be prescribed by a physician and performed by an orthopaedic technician. Under no circumstances should the adjustment made by the physician/orthopaedic technician be altered. Use of the product is recommended for a single patient. With hypersensitive patients, direct contact with skin may cause redness or irritation. In the case of onset of pain, swelling, lumps or any unusual reaction, immediately contact your physician.

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## ENGLISH

**CHARACTERISTICS AND MATERIALS REF. 1090/1091/1093**

- Structure in Sensitive® fabric, elastic, breathable, ladder-proof and antistatic
- Elastic straps for variable and personalised support
- Velcro® front fastening
- 4 posterior guides for insertion of rigid reinforced stays (options) to replace the ones provided
- The abdominal section is reinforced by spiral stays fixed in an oblique position
- Individually adaptable to the patient's pathology

**INDICATIONS REF. 1090/1091/1093**

- Lumbago
- Sciatic lumbalgia and lumbar cruralgia
- Spondylarthrosis and lumbar discopathies
- Slight traumas to the lumbosacral column and paravertebral lumbar contractures
- Post-operative recovery
- REF. 1091 only: After-effects of stable compression fractures of lumbar vertebrae
- REF. 1091 only: Treatment of fractures to transverse processes

### CHARACTERISTICS AND MATERIALS REF. 1092

- Padded shoulder straps
- Structure in Sensitive® fabric, elastic, breathable, ladder-proof and antistatic
- Elastic straps for variable and personalised support
- Velcro® closure on the front
- The abdominal section is reinforced by spiral stays fixed in an oblique position
- Dorso-lumbar part supported by adjustable stays, integrated into the belt
- 4 posterior guides for insertion of rigid reinforced stays (options) to replace the ones provided
- Individually adaptable to the patient's pathology

**INDICATIONS REF. 1092**

- Treatment of back pain caused by diffuse spinal osteoarthritis
- Rheumatism
- Post-surgical protection of vertebral fractures treated surgically with instruments and for dorsal herniotomy

### CONTRAINDICATIONS

Currently no known

### PUTTING ON THE APPLIANCE REF. 1090/1091/1093

- 1 Open the corset and temporarily fasten the Velcro® straps on the side (fig. A).
- 2 Put the corset in place with the loop on the upper side.
- 3 Insert your thumb in the loop and fasten well the central Velcro® closure (fig. B). Make sure the closure is centered.
- 4 Tighten the corset by fastening the straps on the front Velcro® band (fig. C). ONLY on models REF. 1091: first tighten the lower straps, then the top ones.
- 5 If the use of the lumbar pad (REF. 1078/P) is prescribed, fix the pad on the inside to the back Velcro® before putting the corset on.

### PUTTING ON THE APPLIANCE REF. 1092

- 1 Open the corset and temporarily fasten the Velcro® straps on the side. Put the corset in place with the loop on the upper side.
- 2 Insert your thumb in the loop and fasten well the central Velcro® closure (fig. D).
- 3 Make sure the closure is centered. If necessary reposition the corset (fig. E).
- 4 Tighten the corset by fastening the lower straps on the front Velcro® band (fig. F) first.
- 5 Then tighten the upper straps and fasten them on the front Velcro® band (fig. G).
- 6 Place the upper tensioning straps over the shoulders, then feed under the arms and cross over at the back. Make sure they go through the loops (fig. H). Then secure them to the front Velcro® band (fig. I).

**NOTE:** to ensure optimal stabilisation, mould the rigid stays provided and slide them into the pockets provided in the back.

### Attaching the rigid pad:

- if necessary, heat and mould the plastic part to fit the shape of the patient's body.
- place the padding on the plastic part (concave side-without any writing-fig.L)
- apply the adhesive Velcro® fasteners centrally on the opposite side to the padding (fig.M)
- attach the pad to the central Velcro® to the inside of the back of the corset

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