Headquarter: ORTHOSERVICE AG
Via Milano 7 · CH-6830 Chiasso (TI) · Switzerland
Tel. 0041 (0) 91 822 00 88 · Fax 0041 (0) 91 822 00 89
info@orthoservice.com · www.orthoservice.com

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ENGLISH

CHARACTERISTICS AND MATERIALS REF. 1090/1091/1093

- Structure in Sensitive® fabric, elastic, breathable, ladder-proof and antistatic
- Elastic straps for variable and personalised support
- Velcro® front fastening
- 4 posterior guides for insertion of rigid reinforced stays (options) to replace the ones provided
- The abdominal section is reinforced by spiral stays fixed in an oblique position
- Individually adaptable to the patient's pathology

INDICATIONS REF. 1090/1091/1093

- Lumbago
- Sciatic lumbalgia and lumbar cruralgia
- Spondylarthrosis and lumbar discopathies
- Slight traumas to the lumbosacral column and paravertebral lumbar contractures
- Post-operative recovery
- REF. 1091 only: After-effects of stable compression fractures of lumbar vertebrae
- REF. 1091 only: Treatment of fractures to transverse processes

CHARACTERISTICS AND MATERIALS REF. 1092

- Padded shoulder straps
- Structure in Sensitive® fabric, elastic, breathable, ladder-proof and antistatic
- Elastic straps for variable and personalised support
- Velcro® closure on the front
- The abdominal section is reinforced by spiral stays fixed in an oblique position
- Dorso-lumbar part supported by adjustable stays, integrated into the belt
- 4 posterior guides for insertion of rigid reinforced stays (options) to replace the ones provided
- Individually adaptable to the patient's pathology

INDICATIONS REF. 1092

- Treatment of back pain caused by diffuse spinal osteoarthritis
- Rheumatism
- Post-surgical protection of vertebral fractures treated surgically with instruments and for dorsal herniotomy

CONTRAINDICATIONS

Currently no known

PUTTING ON THE APPLIANCE REF. 1090/1091/1093

- 1 Open the corset and temporarily fasten the Velcro® straps on the side (fig. A).
- 2 Put the corset in place with the loop on the upper side.
- 3 Insert your thumb in the loop and fasten well the central Velcro® closure (fig. B). Make sure the closure is centered.
- 4 Tighten the corset by fastening the straps on the front Velcro® band (fig. C). ONLY on models REF. 1091: first tighten the lower straps, then the top ones.
- 5 If the use of the lumbar pad (REF. 1078/P) is prescribed, fix the pad on the inside to the back Velcro® before putting the corset on.

PUTTING ON THE APPLIANCE REF. 1092

- 1 Open the corset and temporarily fasten the Velcro® straps on the side. Put the corset in place with the loop on the upper side.
- 2 Insert your thumb in the loop and fasten well the central Velcro® closure (fig. D).
- 3 Make sure the closure is centered. If necessary reposition the corset (fig. E).
- 4 Tighten the corset by fastening the lower straps on the front Velcro® band (fig. F) first.
- 5 Then tighten the upper straps and fasten them on the front Velcro® band (fig. G).
- 6 Place the upper tensioning straps over the shoulders, then feed under the arms and cross over at the back. Make sure they go through the loops (fig. H). Then secure them to the front Velcro® band (fig. I).

NOTE: to ensure optimal stabilisation, mould the rigid stays provided and slide them into the pockets provided in the back.

Attaching the rigid pad:

- if necessary, heat and mould the plastic part to fit the shape of the patient's body.
- place the padding on the plastic part (concave side-without any writing-fig.L)
- apply the adhesive Velcro® fasteners centrally on the opposite side to the padding (fig.M)
- attach the pad to the central Velcro® to the inside of the back of the corset

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