stabil-T

REF. P1000
3-Punkt-
Hyperextensionsortese
3-Point hyperextension brace
Orthèse d'hyper-extension
à trois points
3-Pistene hyperekstensio-
ortoosi
Iperestensore a tre punti

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ENGLISH

Please read the instructions carefully and thoroughly

DECLARATION OF CONFORMITY

As the manufacturer, ORTHOSERVICE AG declares with sole responsibility that this product is a Class I medical product and has been manufactured in accordance with the basic requirements of the directives 93/42/EEC and 2007/47/EEC. These operating instructions were created in accordance with the specification of the above-mentioned directives. They aim to ensure the safe and correct use of the medical product.

TRADEMARKS OF MATERIALS

Velcro® is a registered trademark of Velcro Industries B.V.

SELECTION/SIZE

Item: REF.P1000	XS		S		MS		M		ML		LS		L		XLS		XL	
Size																		
Pelvic circumf. cm	50/60	60/75	75/90	75/90	75/90	90/105	90/105	105/115	105/115									
Height cm*	30/37	34/42	34/41	36/46	40/52	36/46	40/52	39/49	42/54									
Colour grey																		

*measured from the proximal side of the corpus sterni to the upper edge of the pubic bone (tuberculum pubicum)

OPTION:

PR1-IP35/PC Connect, Connection plate with cervical collar CerviStable
PR1-IP35/P Sternal half-sized pad for pace-maker wearers
PR1-P1000/IA Hypoallergenic non-toxic protective underarm sleeve

INDICATIONS

- Traumatic fractures [T10-L1/L2]
- Osteoporotic or metastatic-related compression fractures of the thoracic spine and the dorsolumbar transition [T10-L1/L2]
- Osteomalacia with fractures
- Dorsolumbar osteoarthritis
- Substitute for plaster corset

MAINTENANCE

- Do not bleach No chemical cleaning
- Do not iron Do not tumble-dry

Washing instructions: hand wash at max. 30° C with neutral soap. Dry away from heat sources.

Do not dispose of the product or any of its components into the environment.

PRECAUTIONS

It is advised that the pressure exerted by the product does not act on parts of the body suffering from wounds, swelling or lumps. The product should not be tightened too firmly so as to avoid the development of pressure points or compression of the nerves or blood vessels underneath. It is advisable to wear the product over a garment, avoiding direct contact with the skin. If in doubt about how to apply, contact an orthopaedic technician. The device should not be worn in the vicinity of open flames or strong electromagnetic fields.

WARNINGS

The product may only be used in the event of medical indication or when prescribed by a physician and must be applied by an orthopaedic technician, the competent authority figure for the application as well as for information on the safe use, according to individual needs. To ensure the effectiveness, tolerability, and correct operation, it is necessary that the device is put on with the utmost care. Any changes to the structure or adjustment of the device must be prescribed by a physician and performed by an orthopaedic technician. Under no circumstances should the adjustment made by the physician/orthopaedic technician be altered. Use of the product is recommended for a single patient. With hypersensitive patients, direct contact with skin may cause redness or irritation. In the case of onset of pain, swelling, lumps or any unusual reaction, immediately contact your physician. Contact an orthopaedic technician for the replacement of used parts. Do not release the tension of the belt abruptly, in order to prevent sudden collapse of the trunk.

WARRANTY

This warranty is governed by the statutory provisions of the country in which the product was purchased. Where relevant that country's regulations concerning warranties and the relationship between the seller and buyer shall apply. If you believe that the terms of the warranty may be applicable, first contact the retailer where the product was purchased. Orthoservice will repair or replace all or part of the device and its accessories if defects in materials and workmanship arise during the period of six months from the date of purchase.

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ENGLISH

CHARACTERISTICS AND MATERIALS

- A lightweight aluminium alloy frame without any sharp metal edges for maximum comfort.
- The brace can be adjusted at the sides with screws and at the top with a telescopic mechanism (the screws do not need to be removed).
- The angle of the underarm joints can be adjusted using the notched system that corresponds with the angle of the underarm pads (7.5° increments).
- An innovative polymer spring enables the pelvic band to be tilted (it can also be secured in position at 0° and at 18° using grub screws).
- Moveable sternal plate.
- A new strap with three-step adjustable length and a clip to hold any excess strap.
- A plastic lever fastening with a safety lock.
- The padded parts are water-repellent and can be immersed in water: PE foam covered with non-toxic synthetic leather.
- The back plate can be worn horizontally and vertically.

ADAPTING FOR THE DOCTOR/TECHNICIAN

- 1 Place the hyperextension brace on the front of the patient's chest (preferably with the patient lying down) to check which adjustments are required (fig. A).
- 2 Loosen the screws on the vertical side stays (fig. B1) and on the horizontal bar (fig. B2).
- 3 Adjust the horizontal bar as follows:
 - Adjust the brace to the width of the patient's chest by moving the two overlapping segments (fig. C1);
 - Tighten the screws, checking that the sternal pad is symmetrical with the sagittal plane (fig. C2);
- 4 Next adjust the vertical stays as follows:
 - Lengthen/shorten the stays so that the pelvic band is level with the pubis and the sternal plate is level with the sternum (fig. D1); remove the screws if necessary, and then replace them.
 - Tighten the screws (fig. D2) and check that the patient can sit down without any discomfort; adjust the vertical stays again if necessary.
 - Check that the two stays are the same length.
- 5 Adjust the angle of the underarm joints (between the vertical stays and the horizontal bar):
 - Remove the screw and the metal cover (fig. E1).
 - Lift the underarm part, turn it to the desired angle, and then reposition it on the notched disc (fig. E2).
 - Replace the metal cover and screws; tighten the screws.
 - Repeat the same procedure on the other side, ensuring that it is adjusted to the same angle.
- 6 Adapt the brace to fit the patient's body exactly:
 - Shape the pelvic band to fit (fig. F1).
 - Shape the underarm pads (fig. F2), ensuring that they fit snugly against the chest area.
 - Shape the side stays at the height of the pelvis (fig. F3), using the crease provided.
- 7 If specified by the doctor, place the lumbar plate against the spine in the position prescribed:
 - Remove the mushroom head screw and position it at the required height (in one of the brass threaded bushings provided) and tighten (fig. G1).
 - Also move the lever to the desired height using the screws provided (fig. G2).
- 8 If specified by the doctor, the lumbar plate can be positioned vertically (fig. H) by removing the strap from the horizontal keepers and passing it through the vertical keepers.
- 9 Adjust the strap as follows:
 - Open the lever fastening, pass the strap behind the patient's back and insert the hook into the mushroom head screw (fig. I1).
 - Adjust the length of the strap so that it provides the right degree of pressure in the lumbar area: slide the strap through the three-step buckle (fig. I2) and fasten by pulling the end (fig. I3); optimum pressure is achieved when the hyperextension brace stays in the correct position even when the lever fastening is open.
 - If the strap is too long, cut it to size or fasten any excess using the clip provided (fig. I4).
 - Click the fastening lever into position (fig. I5).
 - Check that the lumbar plate is centred securely on the spine; fasten it securely to the strap using the keepers provided (fig. I6).

SECURING THE PELVIC BAND (for the doctor/orthopaedic technician)

- The brace is supplied with a pelvic band that tilts and moves. However, it can be secured using the joints at the side if necessary:
- Move the pelvic band so that one of the two holes in the mechanism matches up with the corresponding threaded bushing below (fig. L) until it is at an angle of 0° or 18°.
 - Insert the grub screw provided and tighten (fig. M)
 - Repeat the same procedure on the opposite side.
 - Check that both mechanisms are secured in the same way.

PUTTING ON THE APPLIANCE FOR THE PATIENT

- 1 Open the lever fastening and unhook the strap from the mushroom head screw.
- 2 Place the hyperextension brace on the chest so that the sternal plate rests on the centre of the chest and the pelvic band rests on the pelvis (fig. N)
- 3 Pass the strap behind the back (fig. O) and hook it to the mushroom head screw (fig. P).
- 4 Secure the hyperextension brace using the fastening lever (fig. Q).

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