REF. 1623C *w* 1623H *w*
 1655C *m* 1655H *m*
 1323 *w* 1355 *m*
 1322 *w* 1322SH *w*
 1354 *m* 1622C *w*
 1622H *w* 1654C *m*
 1654H *m* 1321 *w*
 1353 *w*

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CHARACTERISTICS AND MATERIALS

REF.1322-1322SH-1323-1354-1355-1622-1623-1654-1655

- Traction system with elastic straps and loops, for targeted pressure and perfect stabilisation
- 4 posterior guides for insertion of rigid reinforced stays (optional) to replace the ones provided
- Shoulder pads with axillary padding
- Velcro® front fastening
- REF. 1622 - 1623 - 1654 - 1655:
- Posterior insert of reticular fabric for maximum permeability
- REF. 1622 - 1623 - 1654 - 1655: Main fabric made of cotton with a small amount of elastane
- REF. 1322 - 1322SH - 1323 - 1354 - 1355: Fabric made of cotton with a small amount of elastane
- Individually adaptable to the patient's pathology

REF. 1321 - 1353

- Lateral Velcro® closure with accompanying pocket
- Special shaping to favour comfort and adherence
- Ladder-proof hems
- 4 posterior guides for insertion of rigid reinforced stays (optional) to replace the ones provided
- Traction system of wide elastic stretch straps for targeted pressure and perfect stabilisation
- Traction straps with loops, for support and easy wearing
- Velcro® front fastening
- Fabric made of cotton with a small amount of elastane
- Individually adaptable to the patient's pathology

INDICATIONS

- Lumbalgia
- Sciatic lumbalgia and lumbar cruralgia
- Spondylarthrosis and lumbar discopathies
- Slight traumas to the lumbosacral column and paravertebral lumbar contractures
- Post-operation period
- Results of vertebral lumbar relaxations
- Treatment of fracture of transverse processes
- Only REF. 1322 - 1322SH - 1354 - 1622 - 1654 - 1321 - 1353:
- Degenerative lumbar scoliosis
- Spondylosis with or without vertebral listhesis

CONTRAINDIKATIONEN

Currently no known

PUTTING ON THE APPLIANCE


REF.1322-1322SH-1323-1354-1355-1622-1623-1654-1655

- 1 Open up the corset and provisionally affix the Velcro® straps to the strips on the side (fig. A). Position the support in such a way that the inside label faces upwards.
- 2 Apply the support and use the middle Velcro® closure to fix in place securely (fig. B). Ensure that the closure is precisely centred.
- 3 Slowly tighten up the support by securing the lower Velcro® straps to the front strips (fig. C).
- 4 Then secure the upper straps to the front Velcro® strips (fig. D).

REF. 1321 - 1353

- 1 Open up the corset and provisionally affix the Velcro® straps to the strips on the side. Put on the support in such a way that the central Velcro® strip is positioned precisely between the breast bone and the pubic bone (fig. A).
- 2 Guide the side hip strap on the right of the corset over the designated Velcro® closure. To do this, insert fingers into the accompanying pockets (upper loop) and generate the appropriate tension (fig. B).
- 3 Complete securing the side hip strap (lower loop) (fig. C).
- 4 Subsequently tighten the corset and secure the lower tensioning straps, one over the other (left over right, or vice versa) at the height of the pubic bone (fig. D-E).
- 5 Secure the upper tensioning straps in the same way (left over right, or vice versa) and select the required height:
 - a) Overlapping with the lower tensioning straps (fig. F).
 - b) On top of the lower tensioning straps.

To ensure optimal stabilisation, mould the rigid stays provided and slide them into the pockets provided in the back.

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