iper35

REF. IP35

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ENGLISH

Please read the instructions carefully and thoroughly

DECLARATION OF CONFORMITY

As the manufacturer, ORTHOSERVICE AG declares with sole responsibility that this product is a Class I medical product and has been manufactured in accordance with the basic requirements of the directives 93/42/EEC and 2007/47/EEC. These operating instructions were created in accordance with the specification of the above-mentioned directives. They aim to ensure the safe and correct use of the medical product.


TRADEMARKS OF MATERIALS

Velcro® is a registered trademark of Velcro Industries B.V.

SELECTION/SIZE

Size	S	S+	M	M+	L	L+	XL	XL+
Pelvic circumf. cm	60/75	60/75	75/90	75/90	90/105	90/105	105/115	105/115
Height hyperextension brace cm	32/38	38/44	35/41	41/50	37/44	44/52	40/46	46/55

* measured from the proximal side of the corpus sterni to the upper edge of the pubic bone (tuberculum pubicum)
Colour: grey

Dorsolumbal pelotte, can be positioned horizontally and vertically

MAINTENANCE

- Do not bleach
- No chemical cleaning
- Do not iron
- Do not tumble-dry
- Washing instructions: hand wash at max. 30° C with neutral soap (a sponge is recommended); do not dry in the vicinity of heat sources (or, better, dry with a cloth)
- Check screws for tight fit regularly
- Frame: Clean with lukewarm water and neutral soap using a sponge. Dry with a cloth.
- Pads: Remove the pads and clean in lukewarm water with neutral soap. Let dry away from heat sources.
- Do not dispose of the product or any of its components into the environment.
- Contact an orthopaedic technician for the replacement of used parts

PRECAUTIONS

It is advised that the pressure exerted by the product does not act on parts of the body suffering from wounds, swelling or lumps. The product should not be tightened too firmly so as to avoid the development of pressure points or compression of the nerves or blood vessels underneath. It is advisable to wear the product over a garment, avoiding direct contact with the skin. If in doubt about how to apply, contact an orthopaedic technician. The device should not be worn in the vicinity of open flames or strong electromagnetic fields. Do not release the tension of the belt abruptly, in order to prevent sudden collapse of the trunk.

WARNINGS

The product may only be used in the event of medical indication or when prescribed by a physician and must be applied by an orthopaedic technician, the competent authority figure for the application as well as for information on the safe use, according to individual needs. To ensure the effectiveness, tolerability, and correct operation, it is necessary that the device is put on with the utmost care. Any changes to the structure or adjustment of the device must be prescribed by a physician and performed by an orthopaedic technician. Under no circumstances should the adjustment made by the physician/orthopaedic technician be altered. Use of the product is recommended for a single patient. With hypersensitive patients, direct contact with skin may cause redness or irritation. In the case of onset of pain, swelling, lumps or any unusual reaction, immediately contact your physician.

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ENGLISH

CHARACTERISTICS AND MATERIALS

- Telescopic. The screws do not have to be removed, but only loosened for adjustment
- Pelvic band made from flexible polymer. A new fixation mechanism permits a quick and safe change
- Possibility of blocking the spring oscillation using polymer wedges
- Newly designed pin fasteners for easy closing by the patient
- Strap can be adjusted with ratchet closure
- Removable PE foam pads, attached with snap fasteners, covered with synthetic PU material on the outside and blended elastane fabric on the inside
- Base structure made from light aluminium alloy, frame construction
- Movable sternal plate
- Dorsolumbal pelotte, can be positioned horizontally and vertically
- Plus-closure made of light, robust composite material with safety lock
- Protective axilla pads made from hypoallergenic, non-toxic material

ADAPTING FOR THE DOCTOR/TECHNICIAN

- Place the hyperextension orthosis on the front of the patient's torso (if possible the patient should be lying down), to determine which adjustments are to be carried out (fig. A).
- Loosen the screws at the side splint and the front splint (fig. B). Attention: Loosen the screws with max. two turns.
- Proceed to the horizontal adjustment:
 - Adjust the hyperextension orthosis to the thorax of the patient by moving the horizontal splints (fig. C), tighten the screws; check that the sternal pelotte is located symmetrically to the sagittal plane (fig. C);
 - if required, the pelvic band can be exchanged (see separate section).
- Proceed to the vertical adjustment:
 - Shorten/lengthen the hyperextension orthosis until the pelvic band rests on the pubic bone and the plate on the sternum (fig. D); tighten the screws (fig. D); check that the patient can sit down without the orthosis getting in the way. If required repeat the vertical adjustment;
 - check that the side splints are of the same length after the adjustment (this can be effected using the grooves in the metal).
- Carry out precise adjustments to the anatomy of the patient:
 - Model the pelvic band (fig. E). Note: Keep in mind that the ends are made of self-forming polymer;
 - model the axilla parts so that a tight fit at the chest area is ensured (however, they must not exert any pressure) (fig. F); bend the side splints at the height of the orthosis according to the folding mark (fig. G);
 - if required, the retainer plate of the pin fastener can be shaped;
 - if required, the sternal pelotte can also be shaped (fig. H) (in order to achieve a better fit, we recommend using bending irons).
- Attach the back plate on the vertebral column at the place prescribed by the doctor:
 - Loosen the fixing screws of the retainer plate of the pin, move the plate until the desired height is reached and fix again (fig. I);
 - also loosen the fixing screw of the retainer plate of the lever, move until it is at the same height as the pin and fix again (fig. K).
- If prescribed by the doctor, the back plate can also be fixed vertically (fig. L). Pull the strap through the horizontal loops and insert into the vertical loops.
- Adjusting the strap:
 - Ensure that the locking lever is open; loosen the fixing screws of the toothed belt (fig. M) (for easier handling, the bracket can be removed provisionally and only the toothed belt can be used, fig. N)
 - guide the strap over the back of the patient (fig. O); push the toothed belt until the desired length in the ratchet (fig. P); cut any overlapping part of the toothed belt using scissors (fig. P);
 - block the toothed belt by tightening the screws (fig. Q); insert the buckle into the pin (fig. R) and fix the hyperextension orthosis using the locking lever (fig. S)
 - if required, the strap can be readjusted to achieve the desired compressive force in the lumbar region. The optimum compressive force is achieved if the hyperextension orthosis remains in its correct position even when the locking lever is open. Should the strap be too short after the toothed belt has been completely pulled out, it can be extended by cutting the appropriate white Velcro® closures near the lever (fig. T);
 - ensure that the back plate is exactly in the centre of the vertebral column.

INSTRUCTIONS FOR REPLACING THE PELVIC BAND

- Remove the plastic caps at the ends of the pelvic band by removing the respective screws (fig. U and V).
- Pull out the screws at both ends that can be found underneath the caps (fig. W) and remove the pelvic band (fig. X).
- Attach the new pelvic band and fix using the previously removed screws (fig. Y)
- Replace the caps and fix using the appropriate fixing screws (fig. U).

INSTRUCTIONS FOR FIXING THE PELVIC BAND

- If required, the pelvic band, which is supplied movable, can be blocked using the mechanism of the band itself as shown in the following.
- Remove the plastic caps at the ends of the pelvic band by removing the respective screws (fig. U and V).
- Select the wedges with the desired angle (the angle has to be the same on both sides.)
- Move the pelvic band using one hand. With the other hand insert the wedges in the respective receptacles on the mechanism below the caps (fig. Y).
- Replace the caps and fix using the appropriate fixing screws (fig. U).
- Proceed in the same way for the opposite mechanism.

PUTTING ON THE APPLIANCE FOR THE PATIENT

- Open the locking lever and unhook the strap from the pin.
- Place the hyperextension orthosis on the chest (fig. 1), guide the strap over the back of the patient (fig. 2 and 3) and hook into the pin (fig. 4).
- Fit the hyperextension orthosis in such a way that the sternal plate is in the centre of the chest and that the pelvic band does not interfere with sitting down.
- Fix the hyperextension orthosis using the locking lever (fig. 5).

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