GENUFIT 08 / 15  
GENUFIT 08A / 15A

REF. G1108 · G1108A  
REF. G1115 · G1115A  
Kniestütze aus AirX™ Gewebe  
Knee support in AirX™ fabric  
Genouillère en tissu AirX™  
Polvituki AirX-materiaalista  
Ginocchiera in tessuto AirX™

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## ENGLISH

Please read the instructions carefully and thoroughly


### DECLARATION OF CONFORMITY

As the manufacturer, ORTHOSERVICE AG states under its sole responsibility that this product is a Class I medical device that has been manufactured in accordance with the requirements of Directives 93/42/EEC and 2007/47/EC. These instructions have been written in conformity with these directives, and their aim is to ensure the safe and proper use of this medical device.

### SELECTION/SIZE

| Size   | XS    | S     | M     | L     | XL    | XXL   |
|--|-------|-------|-------|-------|-------|-------|
| Thigh circ. taken<br>15 cm above<br>mid-patella cm | 33/39 | 39/45 | 45/52 | 52/58 | 58/66 | 66/75 |
| Length cm  | 34    | 34    | 34    | 34    | 34    | 34    |

Black · Fits right and left

New, very soft and seamless  
popliteal insert made of Soft-X fabric

### MAINTENANCE

- Do not bleach
- No chemical cleaning
- Do not iron
- Do not tumble-dry
- Washing instructions: hand wash at max. 30° C with neutral soap, (REF. G1115 · G1115A only: after removing the hinged bars), do not dry in the vicinity of heat sources

Do not dispose of the product or any of its components into the environment.

### TRADEMARKS OF MATERIALS

Velcro® is a registered trademark of Velcro Industries B.V.

### PRECAUTIONS

It is advised that the pressure exerted by the device does not act on parts of the body suffering from wounds, swelling or lumps. It is advisable to wear a garment, avoiding direct contact with the skin. If in doubt about how to apply, contact an orthopedic technician. The device should not be worn in the vicinity of open flames or strong electromagnetic fields.

### WARNINGS

The product may only be used in the event of medical indication or when prescribed by a physician and must be applied by an orthopaedic technician, the competent authority figure for the application as well as for information on the safe use, according to individual needs. To ensure the effectiveness, tolerability, and correct operation, it is necessary that the device is put on with the utmost care. Any changes to the structure or adjustment of the device must be prescribed by a physician and performed by an orthopaedic technician. Under no circumstances should the adjustment made by the physician/orthopaedic technician be altered. Use of the product is recommended for a single patient. With hypersensitive patients, direct contact with skin may cause redness or irritation. In the case of onset of pain, swelling, lumps or any unusual reaction, immediately contact your physician.

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## ENGLISH

### FEATURES AND MATERIALS

- New anatomical, tight-fit design allows flexible movement and adjusts optimally to the joint
- Breathable, skin-friendly AirX™ fabric
- New, very soft and seamless popliteal insert made of Soft-X fabric for pleasant wearing comfort
- New soft elastic border
- Genufit 08/15 only: Patella opening with foam stabiliser
- Genufit 08A/15A only: New front opening for easier application: internal fastening system for prefixing and external fastening system for final fixing. Patella cut-out
- Redesigned fastening system (European Patent): Stabilizing straps with buckles and elastic insert for the adjustment of optimal pressure settings
- Genufit 08/08A only: New, flexible, flat plastic sleeves for splint and strap attachment
- Genufit 08/08A only: Spiral-shaped lateral support splints
- Genufit 15/15A only: Flat polycentric hinged bars

### INDICATIONS

- REF. G1108 · G1108A only: Initial gonalgia after sprain trauma with mild residual collateral ligament slackness (knee inoperable)
- REF. G1115 · G1115A only: Gonarthrosis or the after-effects of distorsions with residual instability of the collateral ligaments (inoperable knee)
- REF. G1108 · G1115: Mild latero-medial instability of the patella
- REF. G1108 · G1115 only: Femoro-patellar arthrosis

### CONTRAINDICATIONS

Currently no known

### PUTTING ON THE APPLIANCE

- ADAPTING FOR THE DOCTOR/TECHNICIAN**
- 1 Open all Velcro® closures and provisionally close them on themselves.
  - 2 **Genufit 08 · Genufit 15**  
Pull the knee brace over the knee so that the orange label is facing upwards (fig. A) and the patella cut-out is located above the knee cap (fig. B); the guides/joint braces should be parallel to one another.  
**Genufit 08A · Genufit 15A**  
Pull the knee brace on in such a way that the orange label is facing upwards and the patella cut-out is located over the knee cap; then fasten the inner Velcro® closures (fig. C and D) and the outer Velcro® closures (fig. E and F).
  - 3 **Genufit 15 · Genufit 15A**  
If necessary, remove the knee brace, take the joint braces out of teh pockets and bend to correspond to the shape of the leg (fig. G). Then replace the joint braces back in the pockets and put the knee brace back on again.
  - 4 Pull the two Y-shaped ends of the upper strap through the appropriate strap guide, pull fi rmly, simultaneously with both hands and then secure the Velcro® closure (fig. H). Repeat this process for the lower strap (fig. I).

### PUTTING ON THE APPLIANCE FOR THE PATIENT

- 1 Open all Velcro® closures and provisionally close them on themselves.
- 2 **Genufit 08 · Genufit 15**  
Pull the knee brace over the knee so that the orange label is facing upwards (fig. A) and the patella cut-out is located above the knee cap (fig. B); the guides/joint braces should be parallel to one another.  
**Genufi t 08A · Genufi t 15A**  
Pull the knee brace on in such a way that the orange label is facing upwards and the patella cut-out is located over the knee cap; then fasten the inner Velcro® closures (fig. C and D) and the outer Velcro® closures (fig. E and F).
- 3 Pull the two Y-shaped ends of the upper strap through the appropriate strap guide, pull fi rmly, simultaneously with both hands and then secure the Velcro® closure (fig. H). Repeat this process for the lower strap (fig. I).

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