	DE EN FR IT PL DA	Nicht bleichen Do not bleach Ne pas utiliser d'eau de Javel Non candeggiare Nie wybielać Undgå blegemiddel
	DE EN FR IT PL DA	Nicht bügeln Do not iron Ne pas repasser Non stirare Nie prasować Undgå stryging
	DE EN FR IT PL DA	Keine Chemische Reinigung No dry cleaning Pas de nettoyage à sec Non usare detergenti chimici Nie czyścić chemicznie Ingen kemisk rens
	DE EN FR IT PL DA	Nicht im Wäschetrockner trocknen Do not tumble dry Ne pas mettre au sèche-linge Non asciugare in asciugatrice Nie suszyć w suszarce Undgå tørring i tørretumbleren
	DE EN FR IT PL DA	CE-Kennzeichnung Bestätigung der Erfüllung der EU-Anforderungen CE marking Confirmation of compliance with EU requirements Marquage CE Attestation de la conformité aux exigences de l'UE Contrassegno CE Conferma della conformità ai requisiti UE Oznaczenie CE Potwierdzenie spełnienia wymagań UE CE-mærkning Bekræftelse af overholdelse af EU-krav
	DE EN FR IT PL DA	Achtung Lesen Sie alle Warnhinweise und Vorsichtsmaßnahmen in der Gebrauchsanweisung Attention Read all warnings and precautionary measures in the instructions for use Attention Lire toutes les mises en garde et mesures de précaution du mode d'emploi Attenzione Leggere tutte le avvertenze e le misure di precauzione riportate nelle istruzioni per l'uso Uwaga Przeczytaj wszystkie ostrzeżenia i środki ostrożności zawarte w instrukcji obsługi Bemærk Læs alle advarsler og sikkerhedsanvisninger i brugsanvisningen
	DE EN FR IT PL DA	Gebrauchsanweisung Instructions for use Mode d'emploi Istruzioni per l'uso Instrukcja obsługi Brugsanvisning
	DE EN FR IT PL DA	Trocken aufbewahren vor Feuchtigkeit schützen Keep dry protect against moisture Conserver au sec Protéger à l'abri de l'humidité Conservare in un luogo asciutto lontano dall'umidità Przechowywać w suchym miejscu i chronić przed wilgocią Opbevares tørt beskyttes mod fugt
	DE EN FR IT PL DA	Vor Sonnenlicht schützen Protect from sunlight Protéger des rayons du soleil Tenere al riparo dalla luce solare Chronić przed promieniowaniem słonecznym Beskyttes mod sollys
	DE EN FR IT PL DA	Einzelner Patient, mehrfach anwendbar Single patient, multiple use Patient unique, application polyvalente Singolo paziente, uso multiplo Jeden pacjent, wielokrotne użycie En patient, kan anvendes flere gange
	DE EN FR IT PL DA	Webseite zur Patienteninformation Patient information website Site internet d'informations pour les patients Sito web con informazioni per i pazienti Strona internetowa z informacjami dla pacjentów Hjemmeside til patientoplysninger
	DE EN FR IT PL DA	Nach Gebrauch wegwerfen Discard after use Usage unique. Jeter après utilisation Smaltire dopo l'uso Wyrzucić po użyciu Bortskaf efter brug
	DE EN FR IT PL DA	Verwendbar bis Usable until Date d'expiration Utilizzabile fino a Data ważności Kan bruges indtil
	DE EN FR IT PL DA	Hersteller des Produktes Manufacturer of the product Fabricant du produit Produttore dell'articolo Producent produktu Producent
	DE EN FR IT PL DA	Medizinprodukt Medical device Dispositif médical Prodotto medico Urządzenie medyczne Medicinsk produkt
	DE EN FR IT PL DA	Artikelnummer Article number Référence Numero articolo Numer artykułu Varenummer
	DE EN FR IT PL DA	Fertigungslosnummer, Charge Production batch number, batch Numéro de lot de fabrication Numero di lotto di produzione, partita Numer partii produkcyjnej, szarża Batchnummer, batch
	DE EN FR IT PL DA	Pharmazentralnummer Pharmaceutical registration number Numéro central de la pharmacie PZN (Germania) Numer centralny farmacji Farmacentralnummer

SPORLASTIC ORTHOPAEDICS

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INSTRUCTIONS FOR USE EN

INTENDED USE

GENUDYN® OA is a knee brace for relief of the compartment according to the 3-point principle and for physiological guidance and stabilisation.

INDICATIONS

- Conservative/postoperative
- For pain reduction and relief of the medial or lateral compartment in gonarthrosis
- Mild varus/valgus deformities
- All indications and complaints requiring unicompartmental load reduction (e.g. pre- and post-operative for: degenerative meniscus tear, stable tibial plateau fractures, collateral ligament injuries).

CONTRA-INDICATIONS

- Allergic, inflammatory or injury-related skin changes (e.g. swelling, redness) of the body areas to be cared for.
- Impaired circulation or lymphatic soft tissue swelling.
- Neurogenic disorders of the sensory system and skin trophism in the body area to be cared for (sensory disorders with and without skin damage).

SIDE-EFFECTS

If used appropriately and correctly, there are no general side effects, intolerances with long-term effects or allergic reactions known to date.

INTENDED PATIENT TARGET GROUP

The target group is all patients, taking into account indications and contraindications. Instruction and allocation of the correct size is carried out by medical professionals.

PUTTING IT ON AND USE

- **Opening the click locks and releasing the strap tension:** Release the four straps by pressing the respective click lock button on the inside of the brace. Repeat this for all four straps. Press the digit of the Novel button on the outside of the brace and pull on the strap to lengthen it. Repeat this for all four straps **1**.
- **Putting on the brace:** Sit on the edge of a chair and place the brace on the leg with the knee bent at about 30°. The centre of the joint should be at the level of the upper edge of the kneecap (patella). Close the shin strap so that the joint is in the middle of the knee **2**.
- **Adjusting the straps:** Pull the strap 1 tight and make sure that the strap is not twisted. Guide the strap around your lower leg. Close it by inserting it into the click closure on the inside of the brace. A short clicking sound confirms correct closure. Make sure that the strap 1 runs above the calf belly **3**.
- **Tensioning the straps:** You get the necessary strap tension by turning the Novel knob on the outside of the brace. Then repeat these steps for straps 2 - 4. Note that when tightening strap 2, the calf cuff lifts about 1 cm above the shin. This is normal and necessary for a good fit of the brace **4**.
- **Closing the front Velcro fasteners:** The front shin and thigh straps only need to be adjusted by the prosthetist during the first fitting **5**.
- **Check:** Check the position of the brace while standing with the leg extended. The joint centre of the brace should be on the upper edge of the kneecap (patella). You can readjust the straps as needed to ensure a tight fit **6**.
- **Removing the brace:** Press the click lock buttons on the inside of the brace one after the other. You can then pull the brace away to the front **7**.

NOTES FOR THE ORTHOPAEDIC TECHNICIAN

- GENUODYN® OA is delivered with a correct strap sequence for dynamic stabilisation for unicompartamental unloading.
- The fit, joint position and strapping must always be checked for correctness, even when the patient is standing. The front straps remain closed when donning and doffing and the rear straps are opened and closed via the quick-release buckles. All straps should be tightly fastened but should not cause constriction.
- The brace frame is made of aluminium and can be individually bent by the orthopaedic technician.
- All strap lengths are set during fitting and should only be changed later by the patient in case of problems, for example due to circumference changes.

A) NOVEL TURNING WHEEL SYSTEM FOR QUICK ADJUSTMENT:

- The innovative NOVEL system makes it easy to readjust the belt tension - even through clothing **1A**.
- **A** Turning clockwise tightens the belts - the belt tension can be increased in fine doses.
- **B** Press the inner button to release the belt tension.

B) FURTHER ADJUSTMENT AND CUSTOMISATION OPTIONS:

- **Opening the click fasteners:** Release the straps by pressing the click fastener on the inside of the brace. Do this for all four straps **1B**.
- **Releasing the belt tension:** Press the digit on the Novel button and pull on the belt to lengthen it. Do this for all four straps. Next, sit on the edge of a chair and place the brace on the leg at about 30° flexion. The centre of the joint should be at the level of the upper edge of the kneecap (patella). Close the shin strap so that the joint is in the middle of the knee **2B**.
- **Front Velcro fasteners:** The anterior tibia and femur straps do not have click fasteners but a Velcro fastener. Loosen both straps only during the initial fitting. The patient should not readjust these straps **3B**.
- **Closing strap 1:** Pull strap 1 tight and make sure that the strap is not twisted. Guide the strap behind the lower leg and close it with the click closure on the medial side of the brace. Make sure that the strap 1 runs above the calf belly **4B**.
- **Belts 2, 3 and 4:** Repeat the steps from belt 1 **5B**.
- **Close the front Velcro fasteners:** Close the straps so that they fit comfortably and conform to the body and stabilise well **6B**.
- **Check:** Check the position of the brace while standing and with your leg extended. The joint centre of the brace should be on the upper edge of the kneecap (patella) **7B**.
- **Note:** All straps except the tibia and femur strap are numbered consecutively. Make sure that you always close the straps in the correct order. Normally, the tibia and femur straps do not need to be readjusted after the initial fitting **8B**.

C) ADJUSTMENT OF THE CUSHIONS AND BELTS:

- **Pad length:** For optimal fit and wearing comfort, the strap pads must be centred on the leg when the straps are tightened. If a strap pad is too long, remove it from the strap and cut it to fit **1C** with scissors.
- **Strap length:** To adjust the strap length, open the quick release. Loosen the Velcro on the other end of the belt and adjust the length accordingly **2C**.
- **Condylar pads:** The medial and lateral condylar pads are interchangeable and should be adjusted according to the indication. For example, when relieving the medial compartment, the condylar pad should fit the knee laterally and have some free space medially.

D) SETTING THE JOINT LIMIT:

- When delivered, the brace is already set to 0° extension.
- Additional screws and an Allen key for adjusting the extension (EXT) and flexion (FLEX) are included in the scope of delivery.
- **The following joint stops are possible:**
EXT: 20°, 30°, 40°
FLEX: 45°, 60°, 75°, 90°
- 1. **Remove condylar pads:** Remove the condylar pads by detaching them from the inside of the brace (medial and lateral) **1D**.
- 2. **Inserting the screws:** Determine the desired degree of extension and/or flexion and screw the enclosed screws into the corresponding hole using the Allen key provided.
- Screw medial and lateral screws into the same EXT **2D** and/or FLEX **3D** holes. Do not overtighten the screws.
- You can use the side wheel to adjust and fine-tune the corrective force of the brace depending on the degree of gonarthrosis, pain sensation, activity and daily form **4D**.

⚠ IMPORTANT INFORMATION

- Essentially, the indication for and duration of wearing of any orthopaedic device in general, as well as the presence of any of the conditions below in particular, should be discussed with the doctor who is in charge of your treatment.
- Check the functioning of the product with your doctor or specialist retailer.
- The product is designed to treat one patient.
- Do not wear the product directly against open wounds.
- In order for the product to have a long service life and function, it may not be worn in conjunction with oily/greasy or acidic products, ointments or lotions.

- The product can come into contact with splash water.
- Please always close any Velcro fasteners present properly. Damage caused by improper closing of the Velcro fasteners does not constitute a reason for complaint.
- In the event of serious incidents when using this aid, immediately contact your doctor, your health care supply retailer or contact SPORLATIC directly at info@sporlatic.de to report the incident.
- Any improper modifications to the product and/or improper use of the above-mentioned product exclude any product liability on the part of the manufacturer.
- Product may contain metallic objects. We recommend not putting them in contact with electromagnetic interference.
- Possible health-related mutual risks or other disadvantages with certain treatments that may arise in association with the use of the product must be discussed with the doctor in charge of treatment.

 Please return the packaging and the product to the local recycling collection system. Please comply with local regulations.

 SPORLATIC GmbH is affiliated with a returns system for packaging and therefore fully complies with the packaging legislation in the Federal Republic of Germany.

CARE

Cleaning the brace

The pads need to be removed and cleaned regularly depending on the level of activity and perspiration. All pads and cushions are held in place by Velcro fasteners and are therefore easy to remove from the frame, hinges and straps.

Wipe the pads or wash the pads by hand in a mild, antibacterial soapy water (30°C), then rinse well and allow to air dry. Do not dry the product in the tumble dryer, on the heater or in the microwave. Reattach the pads and cushions by pressing them firmly against the corresponding points.

Cleaning the joints

The joints of the brace are pre-lubricated. However, if sand, dirt or water gets into the hinges, they may need to be cleaned and lubricated.

Remove the condylar pad, loosen the hinge screws with the Allen key, clean the hinge thoroughly and then lubricate it with a commercially available synthetic lubricant. Now check that the hinge works properly. Then refit the condylar pad and the hinge screws.

TRANSPORT AND STORAGE CONDITIONS

Please ensure that the product is stored in a dry place and is protected against moisture and sunlight. Store the product at a normal temperature and humidity.

QUALITY MANAGEMENT SYSTEM

All products from SPORLATIC GmbH are subject to product inspection as part of our quality management system. If you still have any complaints about our product, we kindly ask that you contact your specialist retailer.

An online version of the instructions for use can be found on our website at www.sporlatic.de.

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