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 **PLEASE READ THESE INSTRUCTIONS CAREFULLY AND KEEP THEM**

Elbow Brace with R.O.M system

DECLARATION OF CONFORMITY

The manufacturer, ORTHOSERVICE AG declares, under its sole responsibility, that this is a class I medical device that has been manufactured and assessed, according to the EU Regulation 2017/745 (MDR). These instructions have been drawn up in application of the aforementioned Regulation. They are intended to ensure an adequate and safe use of the medical device.

TRADEMARKS OF MATERIALS

Velcro® is a registered trademark of Velcro Industries B.V.

SAFETY PRECAUTIONS

We recommend that the pressure exerted by the device does not act on parts of the body with wounds, swelling, or welts. It is advisable to not over tighten the device so as to avoid generating areas of excessive local pressure or the compression of underlying nerves and/or blood vessels. Furthermore, we recommend wearing a garment to avoid direct contact with the skin. If in doubt on how to apply the device, contact a doctor, physiotherapist, or an orthopedic technician. Read the product composition on the internal label carefully. We recommend to not wear the device in the proximity of free flames or strong electromagnetic fields. Do not apply in direct contact with open wounds.

WARNINGS

It is advisable that the device, designed for the specific indications as set out below, be prescribed by a doctor or a physiotherapist and applied by an orthopedic technician, according to the specific needs of the patient. To ensure effectiveness, tolerability, and proper functioning, the product must be carried out with the utmost care. Never alter the adjustment made by the doctor/physiotherapist/ or orthopedic technician. The responsibilities of the manufacturer shall lapse in case of inappropriate use or adaptation. The orthesis is made and to be used by only one patient. Should the orthesis be inappropriately used, the manufacturer declines all responsibility, as provided by the regulation for medical devices. In hypersensitive individuals, direct contact with the skin may cause redness and irritation. If in doubt on how to apply the device, contact your doctor immediately, and in case of a particularly severe adverse event, report the incident to the manufacturer and the competent authority of your country. The effectiveness of the orthopedic product is only guaranteed when all the components are properly in use.

SELECTION/SIZE

Item	REF.B1014 R / REF.B1014 L	
	S	L
Size		
Arm circ. max. cm	20/28	28/40
Forearm circ. max. cm	20/27	25/35

Color: black-white - Specify right or left

Option: REF. B1015 - Neutral hand support system for Elbo

Option: REF. B1016 - Prono-supination hand support for Elbo

MAINTENANCE

- Do not bleach No chemical cleaning
 Do not iron Do not tumble-dry
 Washing instructions: Do not tumble-dry

- Remove the screws and slide the stay out of the upper arm and forearm cuffs.
- Remove the plastic from the padded parts.
- Hand wash the padded parts in lukewarm water with mild soap; rinse well. Leave to dry away from direct heat.
- Clean the plastic parts with a sponge dipped in lukewarm water and mild soap. Dry with a cloth.

Wash instructions: Do not tumble-dry
Remove the screws and slide the stay out of the upper arm and forearm cuffs. Remove the plastic from the padded parts. Hand wash the padded parts in lukewarm water with mild soap; rinse well. Leave to dry away from direct heat. Clean the plastic parts with a sponge dipped in lukewarm water and mild soap. Dry with a cloth.

Do not use the product or any of its components into the environment.

INDICATIONS

- Stable fractures of elbow, distal humerus, radius and proximal ulna
- Post-surgery immobilization due to tendons and ligaments repair
- Use for restraint after conservative and surgical treatments of elbow fractures
- Elbow luxation and subluxation
- Elbow chronic injuries
- Elbow hyperextension
- Post-trauma treatment

CONTRAINDICATIONS

Currently no known

CHARACTERISTICS AND MATERIALS

- Adjustable position of arm and forearm
- PE plastic to fit with the patient. Pad made of double layer MTP
- Anima interna in plastica (PE) adattabile all'anatomia del paziente, MTP ricoveruto dei due côtés
- Articulated element made of light aluminum alloy and steel junction
- Adjustable extension from 0° to 90° (10° step)
- Adjustable flexion from 0° to 90° (10° step)
- Possible joint block from flexion from 0° to 90° every 10°
- Shoulder strap

DONNING THE BRACE

- Undo all the straps and fasten them back on themselves for the moment.
- Loosen (without removing) the screws that attach the hinged stay to the upper arm cuff and the forearm cuff.
- Loosen the washer (G) to unlock the flexion (F) and extension (E) adjustment mechanisms.
- Place the upper arm cuff on the middle of the arm (fig. 1), pass the straps through the rings provided, and fasten with Velcro (first the distal strap and then the proximal strap).
- Place the lower arm cuff on the middle of the forearm (fig. 2), pass the straps through the rings provided, and fasten with Velcro (first the proximal strap and then the distal strap).
- Slide the hinged stay into the upper arm and forearm cuffs until the joint is aligned with the flexion-extension axis of the elbow (use the epicondyles as a reference point). Tighten the screws provided to lock the hinged stay in position (fig. 3-4).

- angle of flexion and extension (fig. 5) as follows:
 - Adjust the extension by matching up the tab (E) with the desired angle of extension. Turn the lower stay together with the tab if necessary.
 - Adjust the flexion using the tab (F) while keeping the extension tab (E) in position.
- Turn the washer (G) so that the wedge (H) slides into the grooves provided and locks in the flexion and extension adjustment mechanisms (tighten at the bottom).
- For more secure fit, we recommend inserting the Allen screw (I) in the hole provided and tightening at the bottom.

- If required, the protective shoulder strap should also be worn (fig. 6):
- Fasten one of the two ends of the strap to the forearm cuff with Velcro.
 - Pass the strap behind the back, over the shoulder and back to the front.
 - Fasten the other end of the strap to the forearm cuff with Velcro (the strap may be shortened if necessary by cutting and fastening).

SUBSEQUENT WEAR

- Undo all the straps and fasten them back on themselves for the moment.
- Place the upper arm cuff on the upper arm, pass the straps through the rings provided, and fasten with Velcro (first the distal strap and then the proximal strap).
- Place the lower arm cuff on the forearm, pass the straps through the rings provided, and fasten with Velcro (first the proximal strap and then the distal strap).

DONNING THE ACCESSORY FOR NEUTRAL FOREARM ROTATION

- If the PR2-B1015 accessory (fig. 7) is required, it can be done as follows:
- Remove the screw with the washer on the lower part of the forearm cuff (not the screws that lock in the hinged stay).
 - Insert the metal part of the accessory into the forearm cuff (fig. 8) so that the stay slides into the hole provided in the plastic and the threaded insert on the plastic slides into the slot on the accessory. Replace the screw with the washer into the threaded insert.
 - Check that the hand support is in the correct position.
 - Once the patient is wearing the brace, slide in the accessory until the handgrip is optimally positioned. Now tighten the screw to that it is locked into position (fig. 9).

- Fasten the Velcro strap provided parallel to the handgrip to keep the hand in the correct position (fig. 10).
- Use the medium length strap with the Velcro strap provided.

DONNING THE ACCESSORY FOR FOREARM PRONATION OR SUPINATION

- If the PR2-B1016 accessory (fig. 11) is required, it can be done as follows:
- Remove the screw with the washer on the lower part of the forearm cuff (not the screws that lock in the hinged stay).
 - Insert the metal part of the accessory into the forearm cuff (fig. 12) so that the stay slides into the hole provided in the plastic and the threaded insert on the plastic slides into the slot on the accessory. Replace the screw with the washer into the threaded insert.
 - Position the hand support with the concave side facing upwards (fig. 13) if the forearm needs to be supine (with the palm of the hand facing upwards); position the hand support with the concave side facing downwards (fig. 14) if the forearm needs to be prone (with the palm of the hand facing downwards).
 - Once the patient is wearing the brace, slide in the accessory until the hand is resting in the optimal position. Now tighten the screw to that it is locked into position (fig. 15).

Fasten the Velcro strap provided perpendicular to the stay to keep the hand in the correct position (fig. 16).

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Velcro® è un marchio di fabbrica registrato di Velcro Industries B.V.