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PLEASE READ THESE INSTRUCTIONS CAREFULLY AND KEEP THEM

Elbow Brace with R.O.M system

DECLARATION OF CONFORMITY
The manufacturer, ORTHOSERVICE AG declares, under its sole responsibility, that this is a class I medical device that
has been manufactured and assessed, according to the EU Regulation 2017/745 (MDR). These instructions have been
drawn up in application of the aforementioned Regulation. They are intended to ensure an adequate and safe use of the
medical device.

TRADEMARKS OF MATERIALS
Velcro® is a registered trademark of Velcro Industries B.V.
SAFETY/PRECAUTIONS

We recommend that the pressure exerted by the device does not act on parts of the body with wounds, swelling, or welts.
It is advisable to not over tighten the device so as to avoid generating areas of excessive local pressure or the compression
of underlying nerves and/or blood vessels. Furthermore, we recommend wearing a garment to avoid direct contact with
the skin.

WARNINGS

It is advisable that the device, designed specifically for the specific indications as set out below, be prescribed by a doctor or a
physiotherapist and applied by an orthopedic technician, according to the specific needs of the patient.

SELECTION/SIZE

Table with columns: Item Size, REF.B1014 R / REF.B1014 L, and values for Arm circ. max. cm and Forearm circ. max. cm.

Color: black-white - Specify right or left
Option: REF. B1015 - Neutral hand support system for Elbo
Option: REF. B1016 - Prono-supination hand support for Elbo

MAINTENANCE

- Do not bleach, Do not iron, Do not tumble-dry, No chemical cleaning, etc. Washing instructions: Remove the screws and slide the stay out of the upper arm and forearm cuffs.

DO NOT DISPOSE OF THE PRODUCT OR ANY OF ITS COMPONENTS INTO THE ENVIRONMENT.

INDICATIONS

- Stable fractures of elbow, distal humerus, radius and proximal ulna
Post-surgery immobilization due to tendons and ligaments repair
Use for retention after conservative and surgical treatments of elbow fractures
Elbow luxation and subluxation
Elbow chronic injuries

CONTRA-INDICATIONS

- Adjustable position of arm and forearm
PE plastic to fit with the patient. Pad made of double layer MTP
Anima interna in plastica (PE) compatibile à l'anatomie du patient, MTP recouvert des deux côtés

DONNING THE BRACE

- Undo all the straps and fasten them back on themselves for the moment.
Loosen (without removing) the screws that attach the hinged stay to the upper arm cuff and the forearm cuff.
Loosen the washer (G) to unlock the flexion (F) and extension (E) adjustment mechanisms.

SUBSEQUENT WEAR

- Undo all the straps and fasten them back on themselves for the moment.
Place the upper arm cuff on the upper arm, pass the straps through the rings provided, and fasten with Velcro (first the
distal strap and then the proximal strap).

DONNING THE ACCESSORY FOR NEUTRAL FOREARM ROTATION

- Remove the screw with the washer on the lower part of the forearm cuff (not the screws that lock in the hinged stay).
Insert the metal part of the accessory into the forearm cuff (fig. 8) so that the stay slides into the hole provided in the
plastic and the threaded insert on the plastic slides into the slot on the accessory. Replace the screw with the washer into
the threaded insert.

DONNING THE ACCESSORY FOR FOREARM PRONATION OR SUPINATION

- Remove the screw with the washer on the lower part of the forearm cuff (not the screws that lock in the hinged stay).
Insert the metal part of the accessory into the forearm cuff (fig. 12) so that the stay slides into the hole provided in the
plastic and the threaded insert on the plastic slides into the slot on the accessory. Replace the screw with the washer into
the threaded insert.

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