

ERSTE ANWENDUNG
FIRST APPLICATION
PRWEĆIE APPLICATION
PERWIE PRIMOKHENIE
STOSOWANIE PO RAZ
PIERWSZY
PUKEMISHJEET
ELLE PUKEMSKERALLE
HET PRODUCT
ERSTE GEBRUIK
PRIMA ORTEZIS ULIKŠANA
POSTAVLJANJE ORTEZE
(PRVA PRIMJENA)
PRIMA APPLICAZIONE

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FRANZÖSISCH PLEASE READ THESE INSTRUCTIONS CAREFULLY AND KEEP THEM

Brace for hip and coxalgia

DECLARATION OF CONFORMITY
As the manufacturer, ORTHOSERVICE AG declares, under its sole responsibility, that this is a class I medical device that has been manufactured and assessed, according to the EU Regulation UE 2017/745 (MDR). These instructions have been drawn up in application of the aforementioned Regulation. They are intended to ensure an adequate and safe use of the medical device.

TRADEMARKS OF MATERIALS
Velcro® is a registered trademark of Velcro Industries B.V.
SAFETY PRECAUTIONS
We recommend that the pressure exerted by the device does not act on parts of the body with wounds, swelling, or welts. It is advisable to not over tighten the device so as to avoid generating areas of excessive local pressure or the compression of underlying nerves and/or blood vessels. Furthermore, we recommend wearing a garment to avoid contact with the skin. In the case of doubts as to the modes of application, contact a doctor, physiotherapist, or an orthopedic technician. Read the product composition on the internal label carefully. We recommend to not wear the device in the proximity of free flames or direct electromagnetic fields, or of strong electromagnetic fields. Do not apply in direct contact with open wounds.

WARNINGS
It is advisable that the device, designed for the specific indications as set out below, be prescribed by a doctor or a physiotherapist and applied by an orthopedic technician, according to the specific needs of the patient. To ensure effectiveness, tolerability, and proper functioning, the application must be carried out with the utmost care. Never alter the adjustments made by the doctor/physiotherapist/ or orthopedic technician. The responsibilities of the manufacturer should be limited to the use of the device as intended. The orthosis is made and to be used by only one patient. Should the orthosis be inappropriately used, the manufacturer declines all responsibility, as provided by the regulation for medical devices. In hypersensitive individuals, direct contact with the skin may cause redness and irritation. In case of pain, swelling, welts, or any other adverse reaction, contact your doctor immediately, and in case of a particularly severe adverse reaction, immediately contact your doctor. The product is intended for home use. The orthopedic effectiveness of the orthosis is only guaranteed when all the components are properly in use.

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