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PLEASE READ THESE INSTRUCTIONS CAREFULLY AND KEEP THEM

Orthosis for gonarthrosis

DECLARATION OF CONFORMIT Y
As the manufacturer, ORTHOSERVICE AG declares, under its sole responsibility, that this is a class I medical device that has been manufactured and assessed, according to the EU Regulation 2017/745 (MDR). These instructions have been drawn up in application of the aforementioned Regulation. They are intended to ensure an adequate and safe use of the medical device.

TRADEMARKS OF MATERIALS
"Velcro" is a registered trademark of Velcro Industries B.V.
SAFETY PRECAUTIONS
We recommend that the pressure exerted by the device does not act on parts of the body with wounds, swelling, or welts. It is advisable to not too tighten the device so as to avoid generating areas of excessive local pressure or the compression of underlying nerves and/or blood vessels.
Use this orthosis only on the knee of a patient who is in contact with a doctor, physiotherapist, or an orthopedic technician.
Read the product composition on the internal label carefully. The device should not be worn in the vicinity of open flames or strong electromagnetic fields. Do not apply in direct contact with open wounds.

WARNINGS
It is advisable that the device, designed for the specific indications as set out below, be prescribed by a doctor or a physiotherapist and applied by an orthopedic technician, according to the specific needs of the patient.
To ensure effectiveness, tolerability, and proper functioning, the application must be carried out with the utmost care. Never alter the adjustment made by the doctor/physiotherapist/ or orthopedic technician. The responsibilities of the manufacturer shall lapse in case of inappropriate use or adaptation. The orthosis is made and to be used by only one patient. Should the orthosis be inappropriately used, the manufacturer declines all responsibility, as provided by applicable laws.
In case of pain, swelling, welts, or any other adverse reaction, contact your doctor immediately, and in case of a particularly severe adverse event, report the incident to the manufacturer and the competent authority of your country. The effectiveness of the orthopedic product is only guaranteed when all the components are properly in use.

ATTENTION! - The product contains small permanent magnets.

SELECTION/SIZE

Item	REF. 20014 right/left	Xs	S	M	L	XL	XXL
Size							
Circumf. middle of the knee cm	30/33	33/36	36/39	39/43	43/48	48/54	
Height cm	40	40	40	40	40	40	
Colour		grey					

RIGHT = for medial gonarthrosis of the right knee (varus right knee) or for lateral gonarthrosis of the left knee (valgus left knee)
LEFT = for medial gonarthrosis of the left knee (varus left knee) or for lateral gonarthrosis of the right knee (valgus right knee)

Do not bleach No chemical cleaning
 Do not iron No tumble-dry
 Washing instructions: Slide out the metal stay
- Handwash in lukewarm water with mild soap, rinse thoroughly
- Leave to dry in neutral air, away from heat sources
- Slide the metal stay back into place

Dispose of this product and its parts responsibly

INDICATIONS
- Moderate or severe unicompartamental gonarthrosis
- Postoperative support for cartilage or of the meniscus

CONTRAINDICATIONS
Currently no known

CHARACTERISTICS AND MATERIALS
- Patented strap system:
- "allegement" unilateral reduction of compression using a three-point system of load
- the system comes into operation when required: decreased pressure when sitting and tension when standing
- easy adjustable tension using winch system
- easy regulation
- Stretchy PlusComfort fabric: cotton inside for a maximum comfort and polyester outside for maximum resistance
- "leaves" Spindel with wind system
- A single-sided stay made of light aluminum alloy with limit extension of (0°, 10°, 20°, 30°, 40° and flexion of 45°, 60°, 75°, 90°
- Fixing buckles with built-in blocking pull for easy tension, magnet to secure pull after using it.
- Technogel condyle pad to dissipate pressure of straps
- Flexible adjuster straps
- Thigh and calf straps with stretchy element for better comfort when moving and movable pad with Velcro®

PUTTING ON THE APPLIANCE - DONNING THE ORTHOSIS FOR THE FIRST TIME
A Undo all the straps (yellow, blue, magenta, green and white), and release the relevant buckle from its attachment.
B Undo the Velcro® fastening on the fabric part.
C Wrap the brace around the knee with the cuff with the blue and magenta tabs at the top. Fasten the brace with the Velcro® fastening provided and centre the opening over the kneecap (fig. 1).
D Fasten the yellow buckle (calf) to the attachment provided (fig.2a) and use the Velcro® to adjust the length of the strap (fig. 2b). If it is too long, cut off the excess part and reattach the Velcro® fastener provided (fig. 2c).
E Adjust the blue buckle (lower diagonal) to the attachment provided (fig.4a) and use the Velcro® to adjust the length of the strap (fig. 4b). If it is too long, cut off the excess part and reattach the Velcro® fastener provided.
F Fasten the magenta buckle (upper diagonal) to the attachment provided (fig.5a) and use the Velcro® to adjust the length of the strap (fig. 5b). If it is too long, cut off the excess part and reattach the Velcro® fastener provided.
H Bend the leg.
I Attach the white buckle (rear) to the attachment provided (fig. 6a). Adjust the length of the strap using the Velcro® in the middle of the strap (yellow, blue, magenta, green and white), and release the relevant buckle from its attachment.
J Undo the Velcro® fastening on the fabric part.
K Wrap the brace around the knee with the cuff with the blue and magenta tabs at the top. Fasten the brace with the Velcro® fastening provided and centre the opening over the kneecap (fig. 10).
L Fasten the strap using the pull tab provided and fasten the yellow buckle (calf) (fig. 11a). Place the MiniMagnet pull tab in its magnetic housing.
M Tighten the strap using the pull tab provided and fasten the blue buckle (thigh) (fig. 12). Place the MiniMagnet pull tab in its magnetic housing.
N Tighten the strap using the pull tab provided and fasten the green buckle (lower diagonal) (fig. 13). Place the MiniMagnet pull tab in its magnetic housing.
O Tighten the strap using the pull tab provided and fasten the magenta buckle (upper diagonal) (fig. 14). Place the MiniMagnet pull tab in its magnetic housing.
P Keep the leg bent, tighten the strap using the pull tab provided and fasten the white buckle (rear) to the attachment provided (fig. 15). Place the MiniMagnet pull tab in its magnetic housing.

UNDOING THE APPLIANCE - REMOVAL OF THE ORTHOSIS
A Undo all the straps (yellow, blue, magenta, green and white), and release the relevant buckle from its attachment.
B Undo the Velcro® fastening on the fabric part.
C Wrap the brace around the knee with the cuff with the blue and magenta tabs at the top. Fasten the brace with the Velcro® fastening provided and centre the opening over the kneecap (fig. 10).
D Fasten the strap using the pull tab provided and fasten the yellow buckle (calf) (fig. 11a). Place the MiniMagnet pull tab in its magnetic housing.
E Tighten the strap using the pull tab provided and fasten the blue buckle (thigh) (fig. 12). Place the MiniMagnet pull tab in its magnetic housing.
F Tighten the strap using the pull tab provided and fasten the green buckle (lower diagonal) (fig. 13). Place the MiniMagnet pull tab in its magnetic housing.
G Tighten the strap using the pull tab provided and fasten the magenta buckle (upper diagonal) (fig. 14). Place the MiniMagnet pull tab in its magnetic housing.
H Bend the leg.
I Keep the leg bent, tighten the strap using the pull tab provided and fasten the white buckle (rear) to the attachment provided (fig. 15). Place the MiniMagnet pull tab in its magnetic housing.

SUBSEQUENT WEAR
CAUTION: do not undo the Velcro® straps (which were adjusted during the initial fitting), instead use only the levers to fasten the brace.
A Undo all the straps (yellow, blue, magenta, green and white), and release the relevant buckle from its attachment.
B Undo the Velcro® fastening on the fabric part.
C Wrap the brace around the knee with the cuff with the blue and magenta tabs at the top. Fasten the brace with the Velcro® fastening provided and centre the opening over the kneecap (fig. 10).
D Fasten the strap using the pull tab provided and fasten the yellow buckle (calf) (fig. 11a). Place the MiniMagnet pull tab in its magnetic housing.
E Tighten the strap using the pull tab provided and fasten the blue buckle (thigh) (fig. 12). Place the MiniMagnet pull tab in its magnetic housing.
F Tighten the strap using the pull tab provided and fasten the green buckle (lower diagonal) (fig. 13). Place the MiniMagnet pull tab in its magnetic housing.
G Tighten the strap using the pull tab provided and fasten the magenta buckle (upper diagonal) (fig. 14). Place the MiniMagnet pull tab in its magnetic housing.
H Bend the leg.
I Keep the leg bent, tighten the strap using the pull tab provided and fasten the white buckle (rear) to the attachment provided (fig. 15). Place the MiniMagnet pull tab in its magnetic housing.

UNDOING THE APPLIANCE - REMOVAL OF THE ORTHOSIS
A Undo all the straps (yellow, blue, magenta, green and white), and release the relevant buckle from its attachment.
B Undo the Velcro® fastening on the fabric part.
C Wrap the brace around the knee with the cuff with the blue and magenta tabs at the top. Fasten the brace with the Velcro® fastening provided and centre the opening over the kneecap (fig. 10).
D Fasten the strap using the pull tab provided and fasten the yellow buckle (calf) (fig. 11a). Place the MiniMagnet pull tab in its magnetic housing.
E Tighten the strap using the pull tab provided and fasten the blue buckle (thigh) (fig. 12). Place the MiniMagnet pull tab in its magnetic housing.
F Tighten the strap using the pull tab provided and fasten the green buckle (lower diagonal) (fig. 13). Place the MiniMagnet pull tab in its magnetic housing.
G Tighten the strap using the pull tab provided and fasten the magenta buckle (upper diagonal) (fig. 14). Place the MiniMagnet pull tab in its magnetic housing.
H Bend the leg.
I Keep the leg bent, tighten the strap using the pull tab provided and fasten the white buckle (rear) to the attachment provided (fig. 15). Place the MiniMagnet pull tab in its magnetic housing.

UNDOING THE APPLIANCE - REMOVAL OF THE ORTHOSIS
A Undo all the straps (yellow, blue, magenta, green and white), and release the relevant buckle from its attachment.
B Undo the Velcro® fastening on the fabric part.
C Wrap the brace around the knee with the cuff with the blue and magenta tabs at the top. Fasten the brace with the Velcro® fastening provided and centre the opening over the kneecap (fig. 10).
D Fasten the strap using the pull tab provided and fasten the yellow buckle (calf) (fig. 11a). Place the MiniMagnet pull tab in its magnetic housing.
E Tighten the strap using the pull tab provided and fasten the blue buckle (thigh) (fig. 12). Place the MiniMagnet pull tab in its magnetic housing.
F Tighten the strap using the pull tab provided and fasten the green buckle (lower diagonal) (fig. 13). Place the MiniMagnet pull tab in its magnetic housing.
G Tighten the strap using the pull tab provided and fasten the magenta buckle (upper diagonal) (fig. 14). Place the MiniMagnet pull tab in its magnetic housing.
H Bend the leg.
I Keep the leg bent, tighten the strap using the pull tab provided and fasten the white buckle (rear) to the attachment provided (fig. 15). Place the MiniMagnet pull tab in its magnetic housing.

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"Velcro" è un marchio di fabbrica registrato di Velcro Industries B.V.