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PLEASE READ THESE INSTRUCTIONS CAREFULLY AND SAVE THEM

Functional shoulder brace

DECLARATION OF CONFORMITY

As the manufacturer, ORTHOSERVICE AG declares, under its sole responsibility, that this is a class I medical device that has been manufactured and assessed, according to the EU Regulation 2017/745 (MDR). These instructions have been drawn up in application of the aforementioned Regulation. They are intended to ensure an adequate and safe use of the medical device.

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Velcro® is a registered trademark of Velcro Industries B.V.

PRECAUTIONS

We recommend that the pressure exerted by the device does not act on parts of the body with wounds, swelling, or welts. It is advisable to not over tighten the device so as to avoid generating areas of excessive local pressure or the compression of underlying nerves and/or blood vessels. If in doubt on how to apply the device, contact a doctor, physiotherapist, or an orthopedic technician. Read the product information on the internal label carefully. We recommend to not wear the device in the proximity of free flames. Do not apply in contact with open wounds.

WARNINGS

It is advisable that the device, designed for the specific indications as set out below, be prescribed by a doctor or a physiotherapist and applied by an orthopedic technician, according to the specific needs of the patient. To ensure effectiveness, tolerability, and proper functioning, the application must be carried out with the utmost care.

Never alter the adjustment made by the doctor/physiotherapist/ or orthopedic technician.

The responsibilities of the manufacturer shall lapse in case of inappropriate use or adaptation. The manufacturer is not to be used by only one patient.

Should the orthosis be inappropriately used, the manufacturer declines all responsibility, as provided by the regulation for medical devices. In hypersensitive individuals, direct contact with the skin may cause redness and irritation. In case of pain, swelling, welts, or any other adverse reaction, contact your doctor immediately, and in case of a particularly severe adverse reaction, report the incident to the manufacturer and the competent authority of your country. The effectiveness of the orthopedic product is only guaranteed when all of the components are properly in use.

SELECTION/SIZE

Item	REF.20726 R/L					
Size	S	M	L	XL	XXL	
Circ. biceps cm	22/26	26/30	30/35	35/40	40/45	
Colour	grey					

Specify right or left

MAINTENANCE

- Do not bleach No chemical cleaning
- Do not iron Do not tumble-dry
- Washing instructions: hand wash at max. 30° C with neutral soap (a sponge is recommended); do not dry in the vicinity of heat sources.

Do not dispose of the product or any of its components into the environment.

INDICATIONS

- Conservative treatment of instability following anterior or posterior dislocation of the glenohumeral joint.
- Protection of the shoulder during at-risk sports or occupational activities while waiting for surgery or following stabilisation.
- Useful during the rehabilitation of shoulders with painful hyperlaxity or functional instability.
- Use three months after surgery to repair rotator cuff when starting work or sport, including when managing massive ruptures of rotator cuff
- Helpful in managing the rehabilitation of scapular dyskinesia.

CONTRAINDICATIONS

Currently no known

CHARACTERISTICS AND MATERIALS

- BreathSkill fabric structure:
 - Soft and breathable
 - 3D fabric
 - External coverage made of a strong mesh fabr
- Crossover elastics for a perfect joint stabilization
- Lateral adjustable belt; partially elastic and underarm padding
- Adjustable arm with elastics and Velcro® fastening system
- Reinforcement pad in viscoelastic material, internally applicable with Velcro®

PUTTING ON THE APPLIANCE

INITIAL FITTING

- Undo the buckle and undo the Velcro® fastening on the arm.
- Slide the arm into the tubular part, keeping the buckle towards the top and the elasticated arm fastening towards the bottom. Pull the brace up the arm until it almost touches the neck and the shoulder is completely covered. (fig.1)
- Secure the elasticated arm fastening with Velcro® so that the tubular part fits snugly but without blocking the circulation. (fig.2)
- Hold the rear fastening strap, pass it under the opposite armpit (fig.3a) and fasten the buckle at the front (fig.3b).
- Adjust the length of the fastening strap at the front (fig.4a) and at the back (fig.4b), so that the padded part is centred under the armpit: cut off the excess part and secure the strap with the Velcro® fastening.
- CAUTION: The strap must be taut to prevent the brace from moving but not so tight that it blocks the circulation. WE RECOMMEND continuing step by step so that you do not risk cutting off too much.
- Tighten and fasten the two elasticated straps as prescribed by a doctor (both straps can be passed through the upper loop or one through the lower loop and one through the upper loop, depending on which movements require greater restriction).
 - Fasten the upper strap to the lower Velcro® on the arm (fig.5a).
 - Crossover and fasten the lower strap to the upper Velcro® situated near the buckle (fig.5b).
- Move the arm a little a few times to check that the brace is very stable. If necessary, adjust the rear strap.

SUBSEQUENT WEAR

- Undo the buckle and undo the Velcro® fastening on the arm.
- Slide the arm into the tubular part, keeping the buckle towards the top and the elasticated arm fastening towards the bottom. Pull the brace up the arm until it almost touches the neck and the shoulder is completely covered. (fig.10)
- Hold the rear fastening strap, pass it under the opposite armpit (fig.11a) and fasten the buckle at the front (fig.11b).
- Saisir la sangle de fixation à l'arrière, la faire passer sous l'aiselle controlatérale (fig.11a) et fermer à l'avant la boucle à tenon (fig.11b).
- Mette en tension et fermer les deux sangles élastiques:
 - fermer la sangle supérieure sur le velcro inférieur du bras (fig.5a)
 - fermer en croisant la sangle inférieure sur le velcro supérieur, placé à proximité de la boucle de fermeture (fig.12b)

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